

PACO13222

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GALNT17

**Uniprot**

Q6IS24

**Synonyms:**

Williams-Beuren syndrome chromosome region 17; WBSCR17; DKFZp434I2216; DKFZp761D2324; GALNT16; GALNT20; GALNTL3 ;

**Immunogen:**

Human WBSCR17.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A