

WDR13 Antibody



PACO13226

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

WDR13

Uniprot

Q9H1Z4

Synonyms:

WD repeat domain 13; WDR13; DKFZp779C2057; FLJ20563; MG21 ;

Immunogen:

Human WDR13.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A