

WDR19 Antibody



PACO13228

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

WDR19

Uniprot

Q8NEZ3

Synonyms:

WD repeat domain 19; WDR19; FLJ23127; KIAA1638; ORF26; PWDMP ;

Immunogen:

Human WDR19.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A