

DCAF4L2 Antibody



PACO13229

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

DCAF4L2

Uniprot

Q8NA75

Synonyms:

WD repeat domain 21C; WDR21C; FLJ35775; MGC133227 ;

Immunogen:

Human WDR21C.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A