

DCAF10 Antibody



PACO13231

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

DCAF10

Uniprot

Q5QP82

Synonyms:

WD repeat domain 32; WDR32; FLJ23201; MGC10765; MGC131781 ;

Immunogen:

Human WDR32.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A