

PACO13232

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

WDR37

Uniprot

Q9Y2I8

Synonyms:

WD repeat domain 37; WDR37; FLJ40354; KIAA0982; RP11-529L18.2 ;

Immunogen:

Human WDR37.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A