

# WIPF2 Antibody



PACO13242

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

WIPF2

**Uniprot**

Q8TF74

**Synonyms:**

WAS/WASL interacting protein family, member 2; WIPF2; WICH; WIRE ;

**Immunogen:**

Human WIPF2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A