

XAB2 Antibody



PACO13258

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

XAB2

Uniprot

Q9HCS7

Synonyms:

XPA binding protein 2; XAB2; DKFZp762C1015; HCNP; HCRN; NTC90; SYF1 ;

Immunogen:

Human XAB2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A