

XIAP Antibody



PACO13260

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

XIAP

Uniprot

P98170

Synonyms:

X-linked inhibitor of apoptosis; XIAP; API3; BIRC4; ILP1; MIHA; XLP2 ;

Immunogen:

Human XIAP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A