

XIRP2 Antibody



PACO13261

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

XIRP2

Uniprot

A4UGR9

Synonyms:

xin actin-binding repeat containing 2; XIRP2; CMYA3; DKFZp451G055; DKFZp779G1747; FLJ39111; FLJ40200; MYOMAXIN ;

Immunogen:

Human XIRP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A