

XYLT1 Antibody



PACO13270

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

XYLT1

Uniprot

Q86Y38

Synonyms:

xylosyltransferase I; XYLT1; XT-I; XT1 ;

Immunogen:

Human XYLT1.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A