

YBX2 Antibody



PACO13274

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

YBX2

Uniprot

Q9Y2T7

Synonyms:

Y box binding protein 2; YBX2; CONTRIN; CSDA3; DBPC; MGC45104; MSY2 ;

Immunogen:

Human YBX2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A