

YIPF4 Antibody



PACO13276

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

YIPF4

Uniprot

Q9BSR8

Synonyms:

Yip1 domain family, member 4; YIPF4; FinGER4; MGC11061; Nbla11189 ;

Immunogen:

Human YIPF4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A