

YIPF5 Antibody



PACO13277

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

YIPF5

Uniprot

Q969M3

Synonyms:

Yip1 domain family, member 5; YIPF5; DKFZp313L2216; FinGER5; SB140; SMAP-5; SMAP5; YIP1A ;

Immunogen:

Human YIPF5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A