

YPEL1 Antibody



PACO13280

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

YPEL1

Uniprot

O60688

Synonyms:

yippee-like 1 (Drosophila); YPEL1; FKSG3; MGC64992 ;

Immunogen:

Human YPEL1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A