

ZADH2 Antibody



PACO13287

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

ZADH2

Uniprot

Q8N4Q0

Synonyms:

zinc binding alcohol dehydrogenase domain containing 2; ZADH2; MGC45594 ;

Immunogen:

Human ZADH2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A