

ZBP1 Antibody



PACO13289

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ZBP1

Uniprot

Q9H171

Synonyms:

Z-DNA binding protein 1; ZBP1; C20orf183; DAI; DLM-1; DLM1 ;

Immunogen:

Human ZBP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A