

# ZBTB37 Antibody



PACO13290

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ZBTB37

**Uniprot**

Q5TC79

**Synonyms:**

zinc finger and BTB domain containing 37; ZBTB37; D430004I08Rik; MGC2629 ;

**Immunogen:**

Human ZBTB37.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A