

ZBTB7B Antibody



PACO13293

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ZBTB7B

Uniprot

O15156

Synonyms:

zinc finger and BTB domain containing 7B; ZBTB7B; DKFZp686G01254; THPOK; ZBTB15; ZFP67; ZNF857B; c-Krox; hcKrox ;

Immunogen:

Human ZBTB7B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A