

ZC3H4 Antibody



PACO13294

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

ZC3H4

Uniprot

Q9UPT8

Synonyms:

zinc finger CCCH-type containing 4; ZC3H4; C19orf7; KIAA1064 ;

Immunogen:

Human ZC3H4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A