

ZCCHC11 Antibody



PACO13296

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZCCHC11

Uniprot

Q5TAX3

Synonyms:

zinc finger, CCHC domain containing 11; ZCCHC11; PAPD3 ;

Immunogen:

Human ZCCHC11.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A