

ZDHHC5 Antibody



PACO13298

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ZDHHC5

Uniprot

Q9C0B5

Synonyms:

zinc finger, DHHC-type containing 5; ZDHHC5; DKFZp586K0524; KIAA1748; ZNF375 ;

Immunogen:

Human ZDHHC5.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A