

# ZHX2 Antibody



PACO13308

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ZHX2

**Uniprot**

Q9Y6X8

**Synonyms:**

zinc fingers and homeoboxes 2; ZHX2; AFR1; KIAA0854; RAF ;

**Immunogen:**

Human ZHX2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A