

# ZKSCAN3 Antibody



PACO13309

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

ZKSCAN3

**Uniprot**

Q9BRR0

**Synonyms:**

zinc finger with KRAB and SCAN domains 3; ZKSCAN3; FLJ33906; KIAA0426; ZF47; ZFP306; ZNF306; ZNF309; ZSCAN13; Zfp47; dJ874C20.1 ;

**Immunogen:**

Human ZKSCAN3.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A