

# ZMPSTE24 Antibody



PACO13311

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

ZMPSTE24

**Uniprot**

O75844

**Synonyms:**

zinc metallopeptidase (STE24 homolog, *S. cerevisiae*); ZMPSTE24; FACE-1; FACE1; FLJ14968; STE24; Ste24p ;

**Immunogen:**

Human ZMPSTE24.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A