

ZMYND10 Antibody



PACO13312

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

ZMYND10

Uniprot

O75800

Synonyms:

zinc finger, MYND-type containing 10; ZMYND10; BLU; FLU ;

Immunogen:

Human ZMYND10.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A