

ZNF146 Antibody



PACO13315

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNF146

Uniprot

Q15072

Synonyms:

zinc finger protein 146; ZNF146; MGC125660; MGC125661; OZF ;

Immunogen:

Human ZNF146.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A