

ZNF266 Antibody



PACO13324

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

ZNF266

Uniprot

Q14584

Synonyms:

zinc finger protein 266; ZNF266; HZF1 ;

Immunogen:

Human ZNF266.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A