

ZNF3 Antibody



PACO13325

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNF3

Uniprot

P17036

Synonyms:

zinc finger protein 3; ZNF3; A8-51; FLJ20216; HF.12; KOX25; PP838; Zfp113 ;

Immunogen:

Human ZNF3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A