

# ZNF395 Antibody



PACO13329

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

ZNF395

**Uniprot**

Q9H8N7

**Synonyms:**

zinc finger protein 395; ZNF395; DKFZp434K1210; HDBP2; PBF; PRF-1; PRF1; Si-1-8-14 ;

**Immunogen:**

Human ZNF395.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A