

ZNF597 Antibody



PACO13335

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNF597

Uniprot

Q96LX8

Synonyms:

zinc finger protein 597; ZNF597; FLJ33071 ;

Immunogen:

Human ZNF597.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A