

ZNF750 Antibody



PACO13341

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNF750

Uniprot

Q32MQ0

Synonyms:

zinc finger protein 750; ZNF750; FLJ13841; MGC125667; MGC125668; Zfp750 ;

Immunogen:

Human ZNF750.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A