

ZNF839 Antibody



PACO13342

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

ZNF839

Uniprot

A8K0R7

Synonyms:

zinc finger protein 839; ZNF839; C14orf131; FLJ11132 ;

Immunogen:

Human ZNF839.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A