

ZNHIT2 Antibody



PACO13344

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ZNHIT2

Uniprot

Q9UHR6

Synonyms:

zinc finger, HIT type 2; ZNHIT2; C11orf5; FON; MGC120285; MGC120286 ;

Immunogen:

Human ZNHIT2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A