

CDKN1A Antibody



PACO13350

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

CDKN1A

Uniprot

P38936

Synonyms:

cyclin-dependent kinase inhibitor 1A (p21, Cip1); CDKN1A; CAP20; CDKN1; CIP1; MDA-6; P21; SDI1; WAF1; p21CIP1 ;

Immunogen:

Human P21.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A