

PACO13352

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

YWHAG

**Uniprot**

P61981

**Synonyms:**

tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, gamma polypeptide; YWHAG;14-3-3gamma ;

**Immunogen:**

Human 14 3 3GAMMA.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A