

YWHAQ Antibody



PACO13353

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

YWHAQ

Uniprot

P27348

Synonyms:

tyrosine 3-monooxygenase/tryptophan 5-monooxygenase activation protein, theta polypeptide; YWHAQ;14-3-3;1C5; HS1 ;

Immunogen:

Human 14-3-3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A