

ARHGEF9 Antibody



PACO13369

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ARHGEF9

Uniprot

O43307

Synonyms:

Cdc42 guanine nucleotide exchange factor (GEF) 9; ARHGEF9; COLLYBISTIN; HPEM-2; KIAA0424; PEM-2; PEM2 ;

Immunogen:

Human ARHGEF9-Specific.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A