

# MKS1 Antibody



PACO13378

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MKS1

**Uniprot**

Q9NXB0

**Synonyms:**

Meckel syndrome, type 1; MKS1; BBS13; FLJ20345; MES; MKS ;

**Immunogen:**

Human BBS13.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A