

# BLOC1S1 Antibody



PACO13383

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BLOC1S1

**Uniprot**

P78537

**Synonyms:**

biogenesis of lysosomal organelles complex-1, subunit 1; BLOC1S1; BLOS1; FLJ39337; FLJ97089; GCN5L1; MGC87455; MICoA; RT14 ;

**Immunogen:**

Human BLOC1S1-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A