

MCPH1 Antibody



PACO13385

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

MCPH1

Uniprot

A6MCV0

Synonyms:

microcephalin 1; MCPH1; BRIT1; FLJ12847; MCT ;

Immunogen:

Human BRIT1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A