

# NCBP2 Antibody



PACO13390

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

NCBP2

**Uniprot**

P52298

**Synonyms:**

nuclear cap binding protein subunit 2, 20kDa; NCBP2; CBC2; CBP20; NIP1; PIG55 ;

**Immunogen:**

Human CBP20.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A