

ITGAD Antibody



PACO13395

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

ITGAD

Uniprot

Q13349

Synonyms:

integrin, alpha D; ITGAD; ADB2; CD11D; FLJ39841 ;

Immunogen:

Human CD11d.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A