

CEP250 Antibody



PACO13401

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

CEP250

Uniprot

Q9BV73

Synonyms:

centrosomal protein 250kDa; CEP250; C-NAP1; CEP2; CNAP1; MGC88542 ;

Immunogen:

Human CEP250, C-NAP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A