

CHST13 Antibody



PACO13403

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

CHST13

Uniprot

Q8NET6

Synonyms:

carbohydrate (chondroitin 4) sulfotransferase 13; CHST13; C4ST3; MGC119278; MGC119279; MGC119281 ;

Immunogen:

Human CHST13-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A