

# SS18L1 Antibody



PACO13413

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SS18L1

**Uniprot**

O75177

**Synonyms:**

synovial sarcoma translocation gene on chromosome 18-like 1; SS18L1; CREST; KIAA0693; LP2261; MGC26711; MGC78386 ;

**Immunogen:**

Human CREST.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A