

# CTH Antibody



PACO13415

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

CTH

**Uniprot**

P32929

**Synonyms:**

cystathionase (cystathionine gamma-lyase); CTH; MGC9471 ; CSE

**Immunogen:**

Human CSE.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A