

# CYP2E1 Antibody



PACO13422

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

CYP2E1

**Uniprot**

P05181

**Synonyms:**

cytochrome P450, family 2, subfamily E, polypeptide 1; CYP2E1; CPE1; CYP2E; P450-J; P450C2E ;

**Immunogen:**

Human CYP2E1-Specific.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A