

# DFFA Antibody



PACO13430

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

DFFA

**Uniprot**

O00273

**Synonyms:**

DNA fragmentation factor, 45kDa, alpha polypeptide; DFFA; DFF-45; DFF1; ICAD ;

**Immunogen:**

Human DFFA, DFF45.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A