

DOCK1 Antibody



PACO13436

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

DOCK1

Uniprot

Q14185

Synonyms:

dedicator of cytokinesis 1; DOCK1; DOCK180; ced5 ;

Immunogen:

Human DOCK1-Specific.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A